

Measuring radiation dose through the detection of radiation-induced acoustic waves

Susannah Hickling^a, Hao Lei^b, Maritza Hobson^a, Issam El Naqa^{a,b}

^aMcGill University, Montreal, Canada

^bUniversity of Michigan, Ann Arbor, USA



CIRMS 25th Annual Meeting

March 27th, 2017

Dosimetry in radiation therapy

- Variety of dosimetry techniques and protocols required to ensure patients receive the correct dose in the desired location
- Development of advanced delivery techniques presents a dosimetry challenge
- Novel fast and accurate dosimetry techniques are required

Overview of x-ray acoustic computed tomography (XACT)



Bowen et al., Phys Med Biol 36(4), 1991
Xiang et al., Med Phys 40(1), 2013

Theory: Thermoacoustic effect

- Thermoacoustic wave equation:

$$\nabla^2 p(\mathbf{r}, t) - \frac{1}{v_s^2} \frac{\partial^2}{\partial t^2} p(\mathbf{r}, t) = -\frac{\beta}{C_p} \frac{\partial}{\partial t} H(\mathbf{r}, t)$$

p = Pressure

v_s = Speed of sound

β = Isobaric expansion coefficient

C_p = Specific heat capacity

H = Heat energy

Γ = Grüneisen coefficient

Theory: Thermoacoustic effect

- Thermoacoustic wave equation:

$$\nabla^2 p(\mathbf{r}, t) - \frac{1}{v_s^2} \frac{\partial^2}{\partial t^2} p(\mathbf{r}, t) = -\frac{\beta}{C_p} \frac{\partial}{\partial t} H(\mathbf{r}, t)$$

- Initial pressure distribution

$$p(\mathbf{r})|_{t=0} = H(\mathbf{r}) \cdot \Gamma(\mathbf{r}) \quad \text{where} \quad \Gamma = \frac{v_s^2 \beta}{C_p}$$

p = Pressure

v_s = Speed of sound

β = Isobaric expansion coefficient

C_p = Specific heat capacity

H = Heat energy

Γ = Grüneisen coefficient

Theory: Thermoacoustic effect

- Thermoacoustic wave equation:

$$\nabla^2 p(\mathbf{r}, t) - \frac{1}{v_s^2} \frac{\partial^2}{\partial t^2} p(\mathbf{r}, t) = -\frac{\beta}{C_p} \frac{\partial}{\partial t} H(\mathbf{r}, t)$$

p = Pressure

v_s = Speed of sound

β = Isobaric expansion coefficient

C_p = Specific heat capacity

H = Heat energy

Γ = Grüneisen coefficient

- Initial pressure distribution

$$p(\mathbf{r})|_{t=0} = H(\mathbf{r}) \cdot \Gamma(\mathbf{r}) \quad \text{where} \quad \Gamma = \frac{v_s^2 \beta}{C_p}$$

- Relating pressure to dose



$$p(\mathbf{r}) = D(\mathbf{r}) \cdot \rho(\mathbf{r}) \cdot \Gamma(\mathbf{r})$$

Project objective

- Previous work has demonstrated feasibility of using XACT as a dosimetry tool through simulations
- Experimental work has been limited to irradiation of metal blocks
- Present work aims to demonstrate ability of XACT to act as a relative dosimetry technique in water

Experimental set-up



Beam's eye view




Radiation parameters

- 10 MV flattening filter free photon beam produced by clinical radiotherapy linac (Varian TrueBeam)
- At 10 cm depth (imaging plane):
 - Dose per pulse = 0.77 mGy
 - Temperature increase per pulse = 185 nK
- Pulse repetition frequency = 180 Hz
- Pulse length = 4 μ s

Typical transducer signal




Typical transducer signal






Imaging parameters




- Signals acquired at 60 angles
- Image reconstructed from sinogram using back projection algorithm
- XACT images are relative dose images



Simple field: Dose comparison




Simple field: Dose comparison




Root mean square error = 12.2%

Resolution test: XACT image


Field diagram




XACT image



Resolution test: Dose comparison




MLC field: Dose comparison




Root mean square error = 13%

Complex field: XACT image


Field diagram




XACT image



Complex field: Dose comparison



Complex field: Dose comparison



Root mean square error = 10.9%

Outlook: XACT as a dosimeter

- Motivation for XACT dosimetry:
 - Real-time
 - Volumetric
 - No beam perturbation
 - Energy independent
 - Potential to combine with ultrasound imaging
- Future potential applications:
 - 3D dose mapping
 - Small field dosimetry
 - *In vivo* dosimetry
 - Absolute dosimetry

Acknowledgements

Dr. Stephen Davis

Mr. Kyle O'Grady

Dr. Tanner Connell

Dr. Jan Seuntjens

McGill Medical Physics Unit staff and students



S.H. acknowledges partial support by the CREATE Medical Physics Research Training Network grant of the Natural Sciences and Engineering Research Council (Grand number: 432290)

Thank you!

