

A Novel Calorimeter Design for Synchrotron Produced X-ray Beams

I. El Gamal, J.Dessureault, M. McEwen




Introduction

- Synchrotron produced photons are:
 - Broad Spectrum 10 – 150 keV range
 - High dose rate – up to 10 Gy/s
 - Highly collimated
 - Polarized
 - Pulsed - on the order of ps
- Accurate absorbed dose measurements in synchrotron beams lag those in other modalities.
- A dedicated absorbed dose standard has the potential for significantly lower uncertainty.
- Will expand the use of synchrotron beams as a dosimetric research tool.


The BMIT- Beamline

Energy Range	25-150 keV
Source	Super conducting 4.3T wiggler
Nominal field size	22 x1.1 cm
Beamline length	55 m


Measurement position




EBT3 Film Beam Profile Measurements



Beam is flat in
horizontal direction,
varies significantly in
vertical direction



Arrows denote
distance axis
for profiles



Current Methods

- Current methods rely on ion chamber determined air-kerma rates
- Chambers are calibrated using sources with:
 - Lower dose rates ~ 2 mGy/s
 - Different energies
 - Larger field sizes

$$K_{air} = MN_K \prod P_{Corrections}$$


Calorimeter Prototype

- Measures the radiation induced temperature rise.
- Primary measurement, relies only on measurements performed in synchrotron.

$$D = c_{core} \Delta T \prod P_{corrections}$$

Goal: Determine absorbed dose to water at depth for 65 - 140 keV monochromatic beams.


Design



- Vacuum required to reduce heat loss
- Conversion from aluminum to water will be required
- Portable design
- Set-up time < 1 hour


Design

- Measure temperature rise of Al core in vacuum, 40 x 0.5 x 3 mm.
- NTC thermistor chip in glass used, 1.4 x 0.3 mm.
- Aerogel supports used to position the core.




Thermal Simulation

- Thermal results show fast thermistor response to temperature rise.
- Thermistor connecting wires a significant thermal link – can't use copper
- A 2.5 % correction for conductive heat losses, less than water calorimeters in kV beams.





Vacuum Vessel

- Overall Monte Carlo calculated scatter and attenuation correction ranges from **5-15 %** depending on energy.
- No non-Aluminum material in beam path.
- If dose at depth is desired overall correction for all energies is **< 2 %** .



Preliminary Measurements



Future Work

- Overall uncertainty $< 2\%$ for the absorbed dose to water is possible.
- Extensive synchrotron beam measurements to investigate influence quantities.
- Validate theoretical conversion to absorbed dose to water.
- Investigate other applications for this design.

THANK YOU

Islam El Gamal • PhD Candidate • islamelgamal@cmail.carleton.ca



Carleton
UNIVERSITY



Ottawa Medical Physics Institute



Canadian Light Source **Centre canadien de rayonnement synchrotron**

Part or all of the research described in this publication was performed at the Canadian Light Source, a national research facility of the University of Saskatchewan, which is supported by the Canada Foundation for Innovation (CFI), the Natural Sciences and Engineering Research Council (NSERC), the National Research Council (NRC), the Canadian Institutes of Health Research (CIHR), the Government of Saskatchewan, and the University of Saskatchewan.

