

Comparison of Triple-to-Double Coincidence Ratio Liquid Scintillation Counting Activity Determinations of ^{60}Co Using Field Programmable Gate Array and List-Mode Acquired Data



Peyton Lalain¹, Denis Bergeron², Brittany Broder², Sean Jollota¹, Larry DeWerd¹

¹Department of Medical Physics, University of Wisconsin-Madison, Madison, WI ²National Institute of Standards and Technology, Gaithersburg, MD
31st Annual Council on Ionizing Radiation Measurements and Standards
April 29th – May 1st, 2024

TDCR Method

- Triple-to-Double Coincidence Ratio Liquid Scintillation Counting (TDCR LSC) = A primary method for the measurement of activity using ratio of triple-coincidence and double-coincidence count rates
- Allows an absolute measurement of the activity through a free-parameter efficiency
- Three-photomultiplier tube (PMT) system is used, with an LS vial containing activity and LS solution placed in the center

- Emitted scintillation light is collected by the PMTs
- Logic circuitry used to record the collected light as a triple coincidence or a double coincidence




Figure 1. Photo of the TDCR LSC system used by the Radioactivity Group at NIST.

- Using various optical filters, the efficiency of light collection is altered to provide a range of TDCR values

- Mathematical relationship between TDCR and activity (A):

$$TDCR = \frac{R_{Triples}}{R_{Doubles}} = \frac{\epsilon_{Triples}}{\epsilon_{Doubles}} \quad A = \frac{R_{Doubles}}{\epsilon_{Doubles}}$$

- $TDCR$, $R_{Triples}$, and $R_{Doubles}$ are experimentally measured, $\epsilon_{Triples}$ and $\epsilon_{Doubles}$ are determined using MICELLE2 models for efficiency as a function of TDCR value

- We studied the effects of different data acquisition methods, ^{60}Co beta spectrum used, and coincidence resolving time (CRT) on the resulting activity using the TDCR method

TDCR Efficiency Model

- MICELLE2 efficiencies are calculated from the energy spectrum and the energy-dependent quench function
- Three PMTs have equal quantum detection efficiencies
- Counting efficiency for triple coincidences:

$$\epsilon_T = \int_0^{E_{max}} S(E) \left(1 - e^{-\frac{E \cdot Q(E)}{3 \cdot M}}\right)^3 dE$$

- Counting efficiency for logical sum of double coincidences:

$$\epsilon_D = \int_0^{E_{max}} S(E) \left(3 \left(1 - e^{-\frac{E \cdot Q(E)}{3 \cdot M}}\right)^2 - 2 \left(1 - e^{-\frac{E \cdot Q(E)}{3 \cdot M}}\right)^3\right) dE$$

- $S(E)$ is the beta spectrum expected to affect the efficiency, $Q(E)$ is energy dependent quench function, and M is the free parameter

Classical vs. BetaShape ^{60}Co Beta Spectra Study

- Explored the effect of using the classical ^{60}Co beta spectrum compared to the ^{60}Co beta spectrum calculated using BetaShape (BS)
- BetaShape: Developed by LNHB to improve decay data for beta and electron capture emissions using atomic screening and exchange corrections

Table 1			
	Activity [Bq]	Uncertainty [Bq]	Comments
Classical Beta Spectrum	1268.5	2.4	Calculated uncertainty using the standard deviation of the mean for measured values from ^{60}Co Source 2 with no filter. Only Type A uncertainties considered, Type B uncertainties assumed to be constant across measurements.
BetaShape Spectrum	1266.5	2.4	

Table 1: Calculated activities of ^{60}Co using FPGA data (^{60}Co Source 2, no filter (Filter 0), CRT of 150 ns)




Figure 2: Calculated activity [Bq] as a function of the measured TDCR value calculated using FPGA acquired data for ^{60}Co Source 1 and a CRT of 150 ns with different filters applied

Coincidence Resolving Time Study




Figure 3: Effect of CRT on double coincidence counting rates (a) and TDCR values (b) done using list-mode data taken with ^{60}Co source 2 and Filter 0. Post-measurement data analysis available with list-mode data to change CRT and compare observed values.

FPGA vs. List-Mode Data Acquisition Study

- FPGA = Field Programmable Gate Array
 - Array of gates you can program, turning software into hardware
 - Does logic in real-time
 - Must make alterations to dead times and CRTs in real time while measurements are taken because they are hardware adjustments
- In the process of replacing FPGA system with Caen 5724 digitizer
- Collects list-mode data
 - Every measurement gets a timestamp and a peak height, all analysis is done after measurement
 - Can make alterations to lower discriminator level, extending dead time, and coincidence resolving time after physical measurement
 - In this work, only effects of changing CRT were explored

Table 2					
	Activity [Bq]				Comments
	List-Mode Data	Uncertainty	FPGA Data	Uncertainty	
Classical Beta Spectrum	1266.5	3.1	1268.5	2.4	Calculated uncertainty using the standard deviation of the mean for measured values from ^{60}Co Source 2 with no filter and a CRT of 150 ns for both data sets. Only Type A uncertainties considered, Type B uncertainties assumed to be constant across measurements.
BetaShape Spectrum	1264.5	3.1	1266.5	2.4	

Table 2: Calculated activities of ^{60}Co using list-mode data and FPGA data, done using both the classical ^{60}Co beta spectrum and the BS ^{60}Co beta spectrum (^{60}Co Source 2 with no filter (Filter 0) and CRT of 150 ns)

Conclusions

- CRT of at least 50 ns needed to catch all true coincidences
- FGPA and list-mode acquired activity determinations are consistent within 0.16%, but list-mode data uncertainty is slightly larger than FPGA data uncertainty
- Improved beta spectrum does not reduce the efficiency-dependence of the calculated activities
 - Still see different activities with differences in efficiency

Acknowledgements

This research was made possible by the customers of the University of Wisconsin Accredited Dosimetry Calibration Laboratory, whose ongoing business funds research done at the UW Medical Radiation Research Center.

References

- [1] Kossert, Karsten & Marganiec-Gałązka, Justyna & Mougeot, Xavier & Nähle, Ole. (2017). Activity determination of ^{60}Co and the importance of its beta spectrum. *Applied Radiation and Isotopes*. 134. 10.1016/j.apradiso.2017.06.015.
- [2] Mougeot, Xavier. (2017). BetaShape: A new code for improved analytical calculations of beta spectra. *EPJ Web of Conferences*. 146. 12015. 10.1051/epjconf/201714612015.
- [3] Denis E. Bergeron, Jeffrey T. Cessna, Ryan P. Fitzgerald, Lizabeth Laureano-Pérez, Leticia Pibida, Brian E. Zimmerman, Primary standardization of ^{212}Pb activity by liquid scintillation counting, *Applied Radiation and Isotopes*, Volume 190, 2022, 110473, ISSN 0969-8043, 10.1016/j.apradiso.2022.110473.
- [4] Kossert, Karsten, et al. "Uncertainty determination for activity measurements by means of the TDCR method and the CIEMAT/NIST efficiency tracing technique." *Metrologia*, vol. 52, no. 3, 22 May 2015, <https://doi.org/10.1088/0026-1394/52/3/s172>.
- [5] Dutsov, Chavdar, et al. "Evaluation of the accidental coincidence counting rates in TDCR counting." *Nuclear Instruments and Methods in Physics Research Section A: Accelerators, Spectrometers, Detectors and Associated Equipment*, vol. 977, 11 Oct. 2020, p. 164292, <https://doi.org/10.1016/j.nima.2020.164292>.