

Simulation of x-ray imaging devices for regulatory evaluation

Andreu Badal

Division of Imaging, Diagnostics, and Software Reliability
Office of Science and Engineering Laboratories
Center for Devices and Radiological Health
U.S. Food and Drug Administration

31st Annual CIRMS Meeting – April 29, 2024
The Universities at Shady Grove (Rockville, MD)



1. Introduction to the FDA/CDRH/OSEL/DIDSR

- The **Food and Drug Administration** (FDA) is responsible for protecting the public health, in the U.S.A., by ensuring the *safety, efficacy, and security* of human and veterinary drugs, biological products, and medical devices; and by ensuring the safety of our nation's food supply, cosmetics, and products that emit radiation.



- FDA's White Oak campus in Silver Spring, Maryland

DISCLAIMER: The mention of commercial products herein is not to be construed as either an actual or implied endorsement of such products by the U.S. Department of Health and Human Services.



- The **Center for Devices and Radiological Health (CDRH)** is the part of the FDA dedicated to regulating medical and radiation emitting devices.
 - ~2'000 employees, ~25'000 device submissions/year, ~20'000 manufacturers
- The **Office of Science and Engineering Laboratories (OSEL)** is the part of CDRH dedicated to regulatory science research.
 - 3000+ pre-market regulatory consults for CDRH every year.
- The **Division of Imaging, Diagnostics and Software Reliability (DIDSR)** is the part of OSEL dedicated to imaging research. Research programs:
 - Medical Imaging and Diagnostics
 - Digital Pathology
 - Medical Extended Reality
 - Artificial Intelligence / Machine Learning

2. Medical device classification and regulatory evaluation

- FDA employs a risk-based classification paradigm for medical devices:
 - Medical devices are classified and regulated according to degree of risk to the public.
 - Both the technology and intended use and Indications for use inform the risk level.
 - **Intended Use:** general purpose of the device or its function.
 - **Indications for use:** disease or conditions the device will diagnose, treat, prevent, cure or mitigate, including a description of the intended patient population. [21 CFR 814.20(b)(3)(i)]



FDA Total Product Lifecycle activities



- **Pre-market:**

- **Class I** (lowest risk): **most exempt** of FDA notification; follow general controls (Good manufacturing practice, Registration, Adverse Event Reporting, etc.).
- **Class II**: general and specific controls; **Premarket notification [510(k)]**
- **Class III** (highest risk or no predicate): General Controls; **Premarket Approval [PMA]**
- **Pre-Submissions [Qsub]** available to request feedback before submission
- **Investigational device exemption [IDE]**
- **De Novo Petitions** to classify devices as class I or II (class III is default for new devices)

- **Post-market:**

- **Compliance, Inspections, Recalls, Surveillance:** adverse event reporting or medical device reports, trends, and signals from a variety of sources



510(k) Pre-market Notification


- Path to market for the majority of image acquisition (hardware) devices and image processing (software) devices (Class II)
 - Software *in* a Medical Device (SiMD)
 - Software *as* a Medical Device (SaMD)
 - AI/ML software is evaluated essentially in the same way as other software devices, with addition of information on network architecture and training/testing datasets.
- 510 (k) requires determination that the new device [subject device] is **substantially equivalent** to a legally marketed device [predicate device].

[21 CFR 807.92(a)(3)]

<https://www.fda.gov/regulatory-information/search-fda-guidance-documents/510k-program-evaluating-substantial-equivalence-premarket-notifications-510k>

- Performance data needed to support a pre-market submission depends on submission type, device technology, and intended use of the device:
 - **Computational methods** have a great potential to accelerate the development and regulatory evaluation of new imaging devices


A Scientific evidence for medical device regulatory decision-making comes from four different types of models.




Morrison et al. "Advancing regulatory science with computational modeling for medical devices at the FDA's office of science and engineering laboratories." *Frontiers in Medicine* (2019) <https://doi.org/10.3389/fmed.2018.00241>

3. MC-GPU: digital replica of x-ray imaging devices


- To promote the use of in silico imaging in device evaluation, in 2009 we developed the first GPU-accelerated Monte Carlo x-ray tracking code:
 - **MC-GPU**: a version of the PENELOPE photon transport routines that can be efficiently executed in NVIDIA Graphics Processing Units (GPU) using CUDA.
 - X-ray transport through voxelized geometry with delta scattering; binary tree.
 - Open-source code available at: <https://github.com/DIDSR/MCGPU>
 - Badal and Badano, **Accelerating Monte Carlo simulations of photon transport in a voxelized geometry using a massively parallel GPU**, Medical Physics 36, p. 4878-4880 (2009)




• Simulated radiograph



• Original anatomy




• CT recon. no metal



• CT recon. with metal

Monte Carlo loop in the GPU




- Simple flow chart repeated for each x-ray track:
 - $\sim 10^{11}$ x-ray tracks per simulation, as in real acquisition
 - >400 million x-rays/sec in a NVIDIA GeForce RTX 4090
- The quality of the simulation results has to be evaluated with solid *Verification, Validation and Uncertainty quantification* studies:
 - FDA guidance: **Assessing the Credibility of Computational Modeling and Simulation in Medical Device Submissions**: <https://www.fda.gov/regulatory-information/search-fda-guidance-documents/assessing-credibility-computational-modeling-and-simulation-medical-device-submissions>


- **MC-GPU-PET:** extension of the GPU-accelerated Monte Carlo code to simulate PET imaging: <https://github.com/DIDSR/MCGPU-PET>
 - Example: ^{18}F -FDG PET head scan simulated in 43.7 s in 1 GPU.

INPUT

- Voxelized positron source
- Voxelized anatomy (Zubal phantom)
- Setup: Scanner, Acq. & Sinogram



PET RECONSTRUCTION




- Ideal and reconstructed emission maps

Example application:

Virtual Imaging Clinical Trials for Regulatory Evaluation (VICTRE)


- The VICTRE project demonstrated that it is feasible to use computational modeling to reproduce an existing clinical trial.
 - Context of use: comparing the detectability of a solid mass and a calcification cluster in full-field digital mammography (FFDM) and digital breast tomosynthesis (DBT).
 - Images of 3000 breast phantoms were simulated with a model of a Siemens Mammomat Inspiration system using MC-GPU.
 - Detectability was estimated as area under the ROC curve with a channelized Hotelling observer with Laguerre-Gauss channels.
 - Badano et al., [Evaluation of Digital Breast Tomosynthesis as Replacement of Full-Field Digital Mammography: An In-Silico Imaging Trial](#), JAMA Network Open 1, 185474–185474 (2018)
 - Badal et al., [Mammography and breast tomosynthesis simulator for virtual clinical trials](#), Computer Physics Communications 264, 107779 (2021)

- We implemented a complete **computational pipeline** for virtual breast phantom creation, image generation and interpretation.
 - All tools are open-source and available at: <https://github.com/DIDSR/VICTRE>



- Badano et al., JAMA Network Open 1, 185474 (2018)

- Example mammography simulations in the virtual trial:
 - 50 μm voxels; 85 μm pixels (a-Selenium direct detector model)
 - 4 breast glandularity classes and sizes




(a) Fatty

(b) Scattered

(c) Heterogeneous

(d) Dense

- Example **mammography** and **tomosynthesis** simulation with MC-GPU:
 - 25 projections (Siemens protocol); 10^{11} x-rays/projection (~ 1 mGy)
 - <10 min/mammogram ($\sim 4 \times 10^8$ x-ray/sec in an NVIDIA Tesla V100 GPU)



- Mammogram (left) and two DBT slices (right) of a fatty breast with 3 masses inserted at two depths (arrows)
- The results of the VICTRE in silico trial supported the real clinical trial results of slightly superior performance of DBT compared to mammo.

Summary



- The FDA is working towards increasing the use of modeling and simulation in regulatory submissions and evaluation of medical devices.
- Computational imaging methods that simulate medical imaging devices with high realism can be a source of valuable evidence for research and development, and regulatory evaluation, of imaging devices.
 - Solid methods to establish the credibility of computational modeling results needed.
- The CDRH/OSEL develops open-source regulatory science tools that can support regulatory submissions of medical devices.
 - **Regulatory Science Tools Catalog:** <https://cdrh-rst.fda.gov/>

Thank you!



U.S. FOOD & DRUG
ADMINISTRATION

- More information on FDA/CDRH **medical device regulations**:
<https://www.fda.gov/training-and-continuing-education/cdrh-learn>
- CDRH/OSEL **Regulatory Science Tools Catalog**:
<https://cdrh-rst.fda.gov/>