



X-Ray

SPENCER MICKUM | April 30TH, 2024



Agenda



Presenter



Spencer Mickum
Principal Scientist
STERIS Applied Sterilization
Technologies


Agenda

1. X-Ray Background
 - Why do we use x-ray?
2. X-Ray at STERIS
 - Model the Process
3. Radiation Simulation
 - Review Results
4. Conclusions
 - Visualize future goals

X-Ray Background



- X-ray starts as an E-beam, where electrons are generated and accelerated to gain energy
- Electrons are generated in a vacuum environment by the source (electron gun), located at the outer wall of the cavity
- The electrons are then drawn away and accelerated at the specified energy
- At the exit of the accelerator, the cylindrical-shaped beam of high-energy electrons is transported through the beam line from the accelerator to the radiation vault
- Electrons hit at the end a converter plate made of Tantalum where conversion into X-ray (photon) takes place and then hits the product




X-Ray Background



Benefits of X-ray include:


- Directional dose delivery, avoiding losses of isometric irradiation techniques
- Efficient dose delivery, reducing irradiation processing time and correlated temperature rise
- Rapid dose delivery limits oxidative degradation of polymers
- Ability to process to tight dose specifications through improved DUR
- Incremental lap-based dose delivery, offering flexible and precise process definition across a wide range of doses




X-Ray Background



untreated



Gamma around 250 kGy



X-Ray around 250 kGy

Benefits of X-ray include:


- Improved penetration ability of photon energy, like gamma
- Fast and efficient targeted processing that facilitates scale from carton to full pallets of product
- Flexibility – ability to mix different products with different dose requirements in the same irradiation cycle
- Ability to process to tight dose specifications through improved DUR
- Reduced material degradation with reduced processing times and reduce of the color modification when compared to Gamma irradiation.

X-Ray at STERIS



Overview of the x-ray irradiation processes at STERIS

- Several Sites:
 - Chonburi, Thailand,
 - Daniken, Switzerland,
 - Kuala Ketil, Malaysia,
 - Libertyville, Illinois Radiation Technology Center, USA,
 - Suzhou, China,
 - Venlo, The Netherlands
 - Email ast_info@steris.com for more information
- Using IBA Rhodotrons and MEVEX Linear Accelerators





X-Ray at STERIS



Overview of the x-ray irradiation processes at STERIS

- Several Sites:
 - Chonburi, Thailand,
 - Daniken, Switzerland,
 - Kuala Ketil, Malaysia,
 - Libertyville, Illinois Radiation Technology Center, USA,
 - Suzhou, China,
 - **Venlo, The Netherlands**
 - Email ast_info@steris.com for more information
- Using IBA Rhodotrons and MEVEX Linear Accelerators



X-Ray at STERIS




Overview of the x-ray irradiation processes at STERIS Venlo



Radiation Simulation



X-ray simulation



• Courtesy of TRAD

Radiation Simulation



X-ray simulation Studies:

Using TRAD's RayXpert Version 1.9

1. Change in the magnitude of operating scan width required to bring a product's treatment to process capable
2. Review the importance for placing dosimetry within a certain area for repeatable dosimetric measurements.

• Courtesy of TRAD

Radiation Simulation



X-ray simulation 1

Experiment with adjusting the scan width:

Example single homogenous product pallet of 0.20 g/cc





Range of Scan Widths Simulated:

330 cm

300 cm

260 cm

210 cm



Radiation Simulation







X-ray simulation 1

Experiment with adjusting the scan width:

Example single homogenous product pallet of 0.20 g/cc

Range of Scan Widths Simulated:

	DUR
330 cm	1.45
300 cm	1.62
260 cm	2.11
210 cm	5.82



Radiation Simulation



X-ray simulation 1


Experiment with adjusting the scan width:

Example single homogenous product pallet of 0.20 g/cc

Range of Scan Widths Simulated:

	DUR
330 cm	1.45
300 cm	1.62
260 cm	2.11
210 cm	5.82

Review a 5% Variation around DUR



Radiation Simulation



X-ray simulation 1



Experiment with adjusting the scan width:

Example single homogenous product pallet of 0.20 g/cc

Range of Scan Widths Simulated:

	DUR
330 cm	1.45
300 cm	1.62
260 cm	2.11
210 cm	5.82

Review of 5% Variation around DUR



Radiation Simulation



X-ray simulation 1

Experiment with adjusting the scan width:


Example single homogenous product pallet of 0.20 g/cc

Range of Scan Widths Simulated:

	DUR
330 cm	1.45
300 cm	1.62
260 cm	2.11
210 cm	5.82


Review of 5% Variation around DUR

DUR vs Scan Width



Range of 8 cm to accurately simulate DUR

DUR vs Scan Width



Radiation Simulation



X-ray simulation 1

Experiment with adjusting the scan width:

Example single homogenous product pallet of 0.20 g/cc

Range of Scan Widths Simulated:

	DUR
330 cm	1.45
300 cm	1.62
260 cm	2.11
210 cm	5.82

Knowledge of the treatment conditions for the product is critical to reproducing the expected uniformity in simulation

Radiation Simulation



X-ray simulation 2

Experiment with adjusting the position of routine monitoring dosimeter.

Preparation to adopt automation

Also allowing for optimizing position and sensitivity of placement.




Radiation Simulation



X-ray simulation 2

Experiment with adjusting the position of routine monitoring dosimeter.

Vertical - Y-Axis $\pm 2\sigma$ Sensitivity Study



Conclusions



X-ray promises some improvements in economics, material response and processing efficiency.

Some parameters need to be well understood to develop appropriate relationship between processing parameters and operational goals.

- highlighted importance of mirroring

Use of radiation simulation provides STERIS a powerful tool to predict sterilization performance and inform upon the importance of standard operating procedures.

- enhanced flexibility in placement of dosimeters.

Thank You



Questions?

Spencer Mickum

Principal Scientist, Radiation Technology, Applied Sterilization Technologies

STERIS

Email: spencer_mickum@steris.com

Web: www.steris-ast.com